FLEX Your Syringe Options

BAYER E R

Information for Ordering FLEX Consumables





FLEX Syringe Features



Focus on patient safety with added system features



Increase operational efficiency with interoperability and design



Decrease departmental costs with syringe savings



Colored Beacon Technology

A clear indication of the presence of fluid in the syringe with an orange colored identifier

MEDRAD® FluiDots Embedded on the Syringe

Offers additional assurance of fluid in the syringe prior to injection

Ultra-Clear Barrel Design

Further enhances the visualization of fluid to help prevent user error

Reflective Colored Plunger

Plunger color is reflected at the syringe's tip unless 5 mL or more of air is present

Enhanced Syringe Design

Self-rotating syringe docking technology simplifies installation

2D Barcodes

System uses barcodes to determine syringe size, expiration date, if the syringe has already been used

FLEX Your Syringe Options

FLEX syringe options can offer you the flexibility to reduce operational costs:

- > Optimize contrast media use
- Reduce environmental impact and cut the cost of waste disposal

Dual Syringe Utilization Benefits:

Clinical research has shown that simultaneous injection of contrast and saline, in user-selected ratios, can¹:

- Provide more uniform attenuation of the right and left ventricles
- Minimize artifacts by achieving proper attenuation levels
- Visualize the right coronary arteries and right ventricles in a single study by achieving more uniform attenuation

Suggested Add Ons



IBP TS

MEDRAD® Imaging Bulk

Package Transfer Set

50 per box



PRES50

MEDRAD® PRES PressureRated Extension Set
50 per box

Dual 150 mL Syringe Kits



FLEXD-150-SPK
Dual syringe kit with 60" low
pressure T-tubing with prime
tube + 2 fill spikes
20 per box



FLEXD-150-SCS

Dual syringe kit with 60" low pressure T-tubing with prime tube + 1 large saline spike and 1 small contrast spike 20 per box



FLEXD-150-QFT

Dual syringe kit with 60" low pressure T-tubing with prime tube + quick fill tube 20 per box

Single 150 mL Syringe Kits



FLEXS-150-SPK
Single syringe with 60" low pressure tubing with prime tube + 1 fill spike

50 per box



FLEXS-150-QFT
Single syringe with 60" low pressure tubing with prime tube + quick fill tube 50 per box

200 mL Options for Contrast Delivery

A 200 mL syringe option is offered in addition to the 150 mL syringe size for protocols that require higher volumes of contrast.

Catalog No.	Quantity
FLEXD-200-SPK	20 per box
FLEXD-200-SCS	20 per box
FLEXD-200-QFT	20 per box
FLEXS-200-SPK	50 per box
FLEXS-200-QFT	50 per box

References: 1. Kerl JM, Ravenel JG, Nguyen SA, et al. Right heart: split-bolus injection of diluted contrast medium for visualization at coronary CT angiography. Radiology. 2008;247:356-364. 2. Jensen, CT. Dual Flow Contrast Injection for Coronary CTA Improves Visualization of the Right Heart. NASCI Abstract # 05-A-69-NASCI.



At Your Service Every Step of the Way

Our reliable Solutions Delivery teams offer customized service solution options based on the unique needs of your department with an emphasis on adding value for the lifetime of your equipment.

- > Comprehensive hardware and software services
- > Education and support for clinical use
- > Dedicated project management organization
- VirtualCARE® Remote Support provides remote capabilities for advanced diagnostics and real-time support

Order FLEX disposables by contacting our Customer Support Staff (phone: 800-633-7231, Email ricustomercare@bayer.com) or by visiting our online ordering portal at myORDERS.bayer.com.

Bayer reserves the right to modify the specifications and features described herein or to discontinue any product or service identified in this publication at any time without prior notice or obligation. Please contact your authorized Bayer representative for the most current information.

Bayer, the Bayer Cross, MEDRAD Stellant, MEDRAD Stellant FLEX, MEDRAD FluiDots, MEDRAD, Stellant, MEDRAD PRES, Stellant FLEX, FluiDots, and VirtualCARE are trademarks owned by and/or registered to Bayer in the U.S. and/or other countries. Other trademarks and company names mentioned herein are properties of their respective owners and are used herein solely for informational purposes. No relationship or endorsement should be inferred or implied.

© 2018-2020 Bayer. This material may not be reproduced, displayed, modified or distributed without the express prior written consent of Bayer.



Bayer HealthCare LLC 100 Bayer Boulevard P.O. Box 915 Whippany, NJ 07981 U.S.A.

Phone: +1-412-767-2400

+1-800-633-7231 Fax: +1-412-767-4120

More information on radiologysolutions.bayer.com



Manufacturer

Bayer Medical Care Inc. 1 Bayer Drive Indianola, PA 15051-0780 U.S.A

Phone: +1-412-767-2400

+1-800-633-7231

Fax: +1-412-767-4120